

STORAGE AND HANDLING



ULTRA-PURE ACTIN

Catalog Number	Mass	Quantity
160101 - 3 mg	9 mg total powder	3 mg actin protein
160101 - 4 mg	12 mg total powder	4 mg actin protein

STORE IN A COOL,
DRY ENVIRONMENT

For research use only.

Shipping: shipped at ambient temperatures

Storage Conditions: store in a cool, dry environment

Form: desiccated powder (1 mg actin protein is supplied as 3 mg powder with extra mass attributed to trehalose, a lyoprotectant)

Source: rabbit skeletal muscle

Molecular Weight: ~43 kDa

Purity: >99% (SDS-PAGE)

Buffer Conditions Upon Reconstitution: 2 mM Tris-HCl, 0.2 mM CaCl₂, 0.2 mM ATP, 1 mM DTT, and 0.25 M Trehalose (pH 8.0)

Shelf Life: check product label for expiration date

Technical Notes

- store desiccated at 4°C, avoid hot and humid storage conditions
- reconstitute only when ready to use, keep on ice
- perform clarifying spin to remove protein aggregates
- freeze in liquid Nitrogen after reconstitution if desired, avoid repeated freeze-thaw cycles
- visit www.PureSoluble.com/protocols for common actin polymerization protocols

Storage and Handling

Store Ultra-Pure Actin desiccated at 4 °C. Reconstitute only when ready to use by resuspending in ice-cold ultrapure water to 3 mg/ml [actin protein]. Note that 1 mg actin protein is supplied as 3 mg powder (extra mass attributed to trehalose, a lyoprotectant), and reconstitution/dilution should be based on the actin protein concentration. Dilute to 0.4 mg/ml [actin protein] with Actin Working Buffer (5 mM Tris-HCl, 0.2 mM CaCl₂, pH = 8.0) and add ATP to 0.2 mM and DTT to 0.5 mM. Incubate on ice for 1 hour and mix occasionally with gentle vortexing. Next, clarify the actin to remove any protein aggregates. Centrifuge the actin protein at 14k rpm (21k x g) for 15 minutes at 4°C in a microcentrifuge and recover the supernatant on ice. Add sodium azide to 0.05% and store at 4°C. The actin protein is now ready for experimental use and stable under these storage conditions for several weeks.

If desired, the reconstituted actin protein can be aliquoted into smaller experimental batches and flash frozen for later use. Working on ice, aliquot the actin protein into experimental batches not less than 5 ul in volume. Flash freeze the experimental aliquots in liquid Nitrogen and store at -80 °C. Note that it is not recommended to dilute the actin protein prior to freezing. Thaw only when ready to use by placing briefly in a 37°C water bath. Once the actin protein is approximately halfway thawed, remove from the water bath and thaw to completion with gentle flicking. Immediately place thawed actin protein on ice and continue to work on ice. Discard any unused portion of the experimental aliquots to avoid repeated freeze/thaw cycles.

Applications

Ultra-Pure Actin will polymerize into filamentous F-actin when supplemented with KCl and MgCl₂, and kept above its critical concentration. Visit www.PureSoluble.com/protocols for common actin polymerization protocols.



PurSolutions, LLC
111 10th Ave S, Suite 110
Nashville, TN 37203

<https://www.PureSoluble.com>
info@PureSoluble.com
+1-540-560-3411